

Pulsed Infrared/red Light Therapy (**PILT<sup>SM</sup>**)  
**IR2<sup>TM</sup> Therapy System**  
(Large and small animals)

**\$1,289.00 USD**



RevitaVet<sup>TM</sup> utilizes highly rated LEDs and advanced programming to bring you the most effective non-invasive treatment available for your animals. Utilizing 60 infrared and 40 visible red LEDs, the RevitaVet<sup>TM</sup> IR2<sup>TM</sup> Therapy System delivers approximately 3100 mw of power at the rate of greater than 1 joules per minute per cm. sq. to more than 14 sq. in. of treatment area. This versatile device offers you to select any of 7 frequencies setting to allow you to tailor the treatment to a variety of conditions such as:

- Arthritis
- Dysplasia
- Infections
- Musculoskeletal Injuries
- Pain
- Sore and Pulled Muscles
- Wounds

By triggering various photobiological responses, including increasing circulation at the cellular level through a release of nitric oxide, we are able to both trigger and accelerate the body's natural defense and repair mechanisms, and decrease recovery time by as much as 50%, and in some cases even more!



Household Pets



Exotic Animals



Horses



Zoo Animals

## RevitaVet™ IR2 Specs and Features

- **6" x 3.5" = 21 in. sq. treatment area per array**  
**Approx. 3.18 avg. joule delivery per minute per cm. sq.**
- **30 infrared (880 nm - 100 mw peak pwr.) LEDs per array**
- **20 visible red (640 nm - 5 mw peak pwr.) LEDs per array**
- **7 user friendly frequency settings (Nogier) at a 50% duty cycle**
- **1550 mw average power per array**
- **2 LED arrays per unit**
- **Minimal heat output**
- **15 minute treatment times**
- **Auto-off function**
- **Easy Touch™ pushbuttons (low pressure touch operation)**
- **Hands-free operation**
- **Rechargeable battery with smartcharger**
- **Comprehensive user manual**
- **Free tech support**
- **One year limited warranty**

Soft shell turtle had been hit by a car or boat (unsure which) had been being treated with topicals and was getting close to release. It was not eating in captivity so to speed up the release, I used the RevitaVet infrared unit on it and after two treatments on consecutive days, released it after the second treatment. The side by side before and after pictures with just one treatment was pretty amazing!!

Hirschfeld Veterinary Consulting, Inc.  
Beth M. Hirschfeld, DVM, CVA, VMO - President Hollywood, FL USA

We had to open a colon on a dog for severe intestinal obstruction and the dog had been straining for some time and passing blood. In spite of the pain meds etc. the dog appeared very uncomfortable with a heavy respiration coming out of anesthesia. Upon use of the unit on both sides of the abdomen the animal relaxed, was not painful and would walk without much straining. We have used it twice daily and feel it really improved the dog's comfort level. In short I feel this equipment is a beneficial adjunct to our treatment of animals in our facility.

B.J.O'Connor DVM  
O'CONNOR PET HOSPITAL

*See more case studies on website*

## ORDER FORM

Product	Price	Quantity	Total
IR2	\$ 1289.00		
		Promotional code: _____	
<b>Subtotal</b>			\$ _____
State Sales Tax: _____ %			\$ _____

### Shipping/Handling within the continental U.S.

Shipping Method	1 unit	Add'l Units (ea)	
Ground	\$ 15.00	\$ 9.00	
3 Day Select	\$ 26.00	\$ 12.00	
2 <sup>nd</sup> Day	\$ 35.00	\$ 15.00	
Next Day Air	\$ 50.00	\$ 19.00	

*Subject to*

*\*\*Prices  
Change\*\**

**CUSTOMER**

**Total Sale:** \_\_\_\_\_ \$ \_\_\_\_\_

### INFORMATION

“Make Checks Payable to: *“Light Therapy Technologies”*”

**Billing Address: (where credit card bill in sent)**

**Delivery Address (If different from billing):**

**Contact:** \_\_\_\_\_

**Contact:** \_\_\_\_\_

**Company:** \_\_\_\_\_

**Company:** \_\_\_\_\_

**Address:** \_\_\_\_\_

**Address:** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**Fax:** \_\_\_\_\_

**Fax:** \_\_\_\_\_

**E-mail:** \_\_\_\_\_

**E-Mail:** \_\_\_\_\_

**Signature:** \_\_\_\_\_

**Card No.** \_\_\_\_\_

**Exp. Date:** \_\_\_ / \_\_\_ / \_\_\_

**CID:** \_\_\_\_\_

**Total: \$** \_\_\_\_\_