



Pulsed Infrared/red Light Therapy (**PILT<sup>SM</sup>**)

## IR2<sup>TM</sup> Therapy System

(Large and small animals)



Utilizing 60 infrared and 40 visible red LEDs, the RevitaVet<sup>TM</sup> IR2<sup>TM</sup> therapy system delivers approximately 3100 mw of power at the rate of greater than 1 joules per minute per square inch to more than 14 sq. in. of treatment area. This versatile device allows you to select any of 7 frequencies setting to allow you to tailor the treatment to a variety of conditions such as:

- Arthritis
- Dysplasia
- Infections
- Musculoskeletal Injuries
- Pain
- Sore and Pulled Muscles
- Wounds

By triggering various photobiological responses, including increasing circulation at the cellular level through a release of nitric oxide, we are able to both trigger and accelerate the body's natural defense and repair mechanisms, and decrease recovery time by as much as 50%, and in some cases even more!



Hocks



Hoofs



Tendon/suspensories



Loins/Back



## RevitaVet™ IR2 Specs and Features

- **6" x 3.5" = 21 in. sq. treatment area per array**  
**Approx. 3.18 avg. joule delivery per minute per cm. sq.**
- **30 infrared (880 nm - 100 mw peak pwr.) LEDs per array**
- **20 visible red (640 nm - 5 mw peak pwr.) LEDs per array**
- **7 user friendly frequency settings (Nogier) at a 50% duty cycle**
- **1550 mw average power per array**
- **2 LED arrays per unit**
- **Minimal heat output**
- **15 minute treatment times**
- **Auto-off function**
- **Easy Touch™ pushbuttons (low pressure touch operation)**
- **Hands-free operation**
- **Rechargeable battery with smartcharger**
- **Comprehensive user manual**
- **Free tech support**
- **One year limited warranty**

"We have been using the Revitavet IR2 for the past 3 months and have seen incredible results. We are now using it pre-race on all our horses before sending them to the paddock and as daily therapeutic therapy. One particular case was a 2yr old thoroughbred in race training that become sore in the stifles and over her back. The filly was traveling lame and not wanting to train, and we contemplated taking her home and turning her out to recover, but decided to give light therapy a try with the Revitavet IR2. I targeted the stifles and acupuncture areas throughout the back and within a few days, she had made a dramatic turn around and looked forward to going to the track. She went on and placed in her maiden debut."

Ross Russell, Thoroughbred trainer and Michigan State University Vet student

Two year Quarter horse filly with non-healing laceration on the left carpus. The carpus was swollen, sticking out 2 ½ inches beyond normal and a non-healing flesh wound 3 ½ " in diameter for a four month duration. I applied the RevitaVet to the lesion four days in a row. After the first application I could see the new skin growth. By the fifth day the flesh wound was only 1 ½ " in diameter. In addition, the swelling had reduced. Two weeks later the flesh wound was less than a ½ " in diameter – Four applications

Twenty year old Quarter horse mare with a fractured medial splint bone and soft tissue damage from the carpus distal to two inches above the fetlock. She had caught the leg in a wooden fence. The skin was not broken. She had grade 3 lameness. We applied the RevitaVet on her leg three times, one week apart. The problem was 3 ½ weeks old when I first saw it. There had been no improvement since the incident. After the first application, the swelling was reduced and the lameness was a grade 2. At the time of the second application, she had a grade 1 lameness with reduced swelling. By the third application, the mare was normal with very little swelling. She was returned to regular riding activity. The response after the first treatment was remarkable in light of the previous 3 ½ weeks with no improvement. The mare is doing fine thanks to RevitaVet.

Glen Wilkinson, DVM - Wilkinson Veterinary Clinic, Fremont, TX



**ORDER FORM**

Product	Price	Quantity	Total
IR2	\$ 1289.00		
Subtotal			\$
AZ State Sales Tax ( if applicable)			\$

**Shipping/Handling within the continental U.S.**

Shipping Method	1 unit	Add'l Units (ea)	
Ground	\$ 15.00	\$ 9.00	
3 Day Select	\$ 29.00	\$ 12.00	
2 <sup>nd</sup> Day	\$ 36.00	\$ 15.00	
Next Day Air	\$ 50.00	\$ 19.00	

**Total Sale:** \$ \_\_\_\_\_

**CUSTOMER INFORMATION**

**Personal or Company Checks Accepted** - Make Checks Payable to: *“Light Therapy Technologies”*

<b>Billing Address: (where credit card is sent)</b>	<b>Delivery Address (If different from billing):</b>
Contact: _____	Contact: _____
Company: _____	Company: _____
Address: _____	Address: _____
_____	_____
Phone: _____	Phone: _____
Fax: _____	Fax: _____
E-mail : _____	

Signature: \_\_\_\_\_

Card No. \_\_\_\_\_

Exp. Date: \_\_\_/\_\_\_/\_\_\_      CID: \_\_\_\_\_      Total: \$ \_\_\_\_\_