

IR2 Infrared Therapy and Rehab System



RevitaVet™ IR2 Specs and Features

- 5.5" x 3.5" = 19.25 sq. in. treatment area per array X2 = 38.5 sq. in
- Approx. 3.18 avg. joule delivery per minute per cm. sq.
- 30 infrared (880 nm - 100 mw peak pwr.) LEDs per array
- 20 visible red (640 nm - 5 mw peak pwr.) LEDs per array
- 7 user friendly frequency settings (Nogier) at a 50% duty cycle
- 1550 mw average power per array
- 2 - LED arrays (treating pads) per unit with minimal heat output
- 15 minute treatment times
- 1 - Controller with Auto-off function
- Easy Touch™ pushbuttons (low pressure touch operation)
- 2 - 1/2 inch elastic Velcro straps to secure treating pads on legs
- 1 - 7 foot elastic Velcro strap to secure treating pad on horse's body
- Hands-free operation
- 1 - Portability - rechargeable battery and smart charger set
- 1 - Battery storage bag/holster while using on horse
- 1 - Comprehensive user manual
- 1 - Carry case
- 1 - AC Power supply
- Free tech support

30 Day Money Back or One year limited warranty